smiths medical

bringing technology to life

LABELLING LAYOUT SHEET

L6397

ISSUE - 100

Title: Instructions For Use - Blue Line Ultra® Tracheostomy Inner Cannula,
Single Patient Use - (004/005/248A) (Form Fill Seal - Soft Pack Only)

01 English	12 French	13 German	14 Italian	15 Spanish	16 Portuguese	17 Dutch	18 Swedish	19 Danish	20 Norwegian
Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
21 Finnish	22 Greek	23 Japanese	24 Czech	25 Slovakian	26 Hungarian	27 Polish	28 Slovenian	29 Estonian	30 Latvian
Y	Y	N	Y	N	Y	Y	N	Y	Y
31 Lithuanian	32 Turkish	33 Bulgarian	34 Romanian	35 Russian	36 Croatian	37 Chinese	38 Korean	39 Arabic	
N	Y	Y	N	N	N	N	N	N	

CCN: 004/005/248A

Size: Internal pages: 210mm x 148mm (A5), Folded to 105mm x 148mm (A6).

External pages: 315mm x 148mm, Folded to 105mm x 148mm (A6).

Print Details: PMS Pantone Process Black U.

Design Authority: Hylke (MHYT).

Material: 80gsm White Bond Paper.

Description: Printed Instruction Booklet, 6 pages external cover and 60 pages internal. (See folding diagram on page 2).

NOTES:

A) The layouts comply with the requirements of DVP802 REV. 000.

B) This artwork is based on CCN. 004/005/192 (L5306 Issue 001)

Last Change: PWR 6097, PCO HYT008598

- Introduction of new 'CZ' codes into manufacturing at the Smiths Medical Czech Republic site, Hranice
- Manufacturer address changed to Czech Republic address.
- Catalogue number and CZ Internal Ref added.
- Trademark statement updated.
- Glossary updated






(The following changes are confirmed by Audra Howard Grigg (15 May 14)

- Inner Cannula I.D. (mm) sizes on Table grid changed from 6.5 to 6.0, 6.0 to 6.5, 6.5 to 7.0, 7.0 to 7.5 and 7.5 to 8.0.

- Changes to Note 2 under the Table Grid and Warnings section number 10.

"The use of an Inner Cannula reduces the airway diameter by 1.5mm (excluding size 6.0mm tracheostomy tube which reduces the airway diameter by 1.0mm)" changed to

"The use of an Inner Cannula reduces the airway diameter by 1.0mm (excluding size 10.0mm tracheostomy tube which reduces the airway diameter by 1.6mm)".

SMI Labelling  Name: Emma Howard Date: 16 May 2014 Signature held on master	Regulatory Affairs  Name: Marianne Slade Date: 26 May 2014 Signature held on master	Marketing  Name:  Date:  Signature held on master
Compiled by Dan Johnson 	Design Authority R&D / Operations Technology  Name:  Date:  Signature held on master	Date Effective 13 JUN 2014
Compilation Date - 16 MAY 14	Circulation: Czech Republic	Sheet 1 of 19

LABELLING LAYOUT SHEET

Title: Instructions For Use - Blue Line Ultra® Tracheostomy Inner Cannula,
Single Patient Use - (004/005/248A) (Form Fill Seal - Soft Pack Only)

en English	fr French	de German	it Italian	es Spanish	pt Portuguese	nl Dutch	sv Swedish	da Danish	no Norwegian
Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
fi Finnish	el Greek	ja Japanese	cs Czech	sk Slovakian	hu Hungarian	pl Polish	sl Slovenian	et Estonian	lv Latvian
Y	Y	N	Y	N	Y	Y	N	Y	Y
lt Lithuanian	tr Turkish	bg Bulgarian	ro Romanian	ru Russian	hr Croatian	zh Chinese	ko Korean	ar Arabic	
N	Y	Y	N	N	N	N	N	N	

CCN: 004/005/248A

Size: Internal pages: 210mm x 148mm (A5), Folded to 105mm x 148mm (A6).

External pages: 315mm x 148mm, Folded to 105mm x 148mm (A6).

Print Details: PMS Pantone Process Black U.

Design Authority: Hythe (MHT).

Material: 80gsm White Bond Paper.

Description: Printed Instruction Booklet, 6 pages external cover and 60 pages internal. (See folding diagram on page 2).

NOTES:

A) The layouts comply with the requirements of DVP802 REV. 000.

B) This artwork is based on CCN. 004/005/192 (L5306 Issue 001)

Last Change; PWR 5097, PCO HYT008596

- Introduction of new 'CZ' codes into manufacturing at the Smiths Medical Czech Republic site, Hranice
- Manufacturer address changed to Czech Republic address.
- Catalogue number and CZ Internal Ref added.
- Trademark statement updated.
- Glossary updated



(The following changes are confirmed by Audra Howard Grigg (15 May 14)

- Inner Cannula I.D. (mm) sizes on Table grid changed from 5.5 to 6.0, 6.0 to 6.5, 6.5 to 7.0, 7.0 to 7.5 and 7.5 to 8.0.

- Changes to Note 2 under the Table Grid and Warnings section number 10.

"The use of an Inner Cannula reduces the airway diameter by 1.5mm (excluding size 6.0mm tracheostomy tube which reduces the airway diameter by 1.0mm)" changed to

"The use of an Inner Cannula reduces the airway diameter by 1.0mm (excluding size 10.0mm tracheostomy tube which reduces the airway diameter by 1.5mm)".

SMI Labelling  Name: Emma Howard Date: 16 May 2014. <i>Signature held on master</i>	Regulatory Affairs  Name: Date: <i>Signature held on master</i>	Marketing  Name: DANIEL JARRARD Date: 16 MAY 2014. <i>Signature held on master</i>
Compiled by Dan Johnson 	Design Authority R&D / Operations Technology  Name: Date: <i>Signature held on master</i>	Date Effective 13 JUN 2014
Compilation Date - 16 MAY 14	Circulation: Czech Republic	Sheet 1 of 19

LABELLING LAYOUT SHEET

L6397

ISSUE - 100

Title: Instructions For Use - Blue Line Ultra® Tracheostomy Inner Cannula,
Single Patient Use - (004/005/248A) (Form Fill Seal - Soft Pack Only)

en English	fr French	de German	it Italian	es Spanish	pt Portuguese	nl Dutch	sv Swedish	da Danish	no Norwegian
Y	Y	Y	Y	Y	Y	Y	Y	Y	Y
fi Finnish	el Greek	ja Japanese	cs Czech	sk Slovakian	hu Hungarian	pl Polish	sl Slovenian	et Estonian	lv Latvian
Y	Y	N	Y	N	Y	Y	N	Y	Y
lt Lithuanian	tr Turkish	bg Bulgarian	ro Romanian	ru Russian	hr Croatian	zh Chinese	ko Korean	ar Arabic	
N	Y	Y	N	N	N	N	N	N	

CCN: 004/005/248A

Size: Internal pages: 210mm x 148mm (A5), Folded to 105mm x 148mm (A6).

External pages: 315mm x 148mm, Folded to 105mm x 148mm (A6).

Print Details: PMS Pantone Process Black U.

Design Authority: Hythe (MHT).

Material: 80gsm White Bond Paper.

Description: Printed Instruction Booklet, 6 pages external cover and 60 pages internal. (See folding diagram on page 2).

NOTES:

A) The layouts comply with the requirements of DVP802 REV. 000.

B) This artwork is based on CCN. 004/005/192 (L5306 Issue 001)

Last Change: PWR 5097, PCO HYT008596

- Introduction of new 'CZ' codes into manufacturing at the Smiths Medical Czech Republic site, Hranice
- Manufacturer address changed to Czech Republic address.
- Catalogue number and CZ Internal Ref added.
- Trademark statement updated.
- Glossary updated




(The following changes are confirmed by Audra Howard Grigg (15 May 14)

- Inner Cannula I.D. (mm) sizes on Table grid changed from 5.5 to 6.0, 6.0 to 6.5, 6.5 to 7.0, 7.0 to 7.5 and 7.5 to 8.0.


- Changes to Note 2 under the Table Grid and Warnings section number 10.

"The use of an Inner Cannula reduces the airway diameter by 1.5mm (excluding size 6.0mm tracheostomy tube which reduces the airway diameter by 1.0mm)" changed to



"The use of an Inner Cannula reduces the airway diameter by 1.0mm (excluding size 10.0mm tracheostomy tube which reduces the airway diameter by 1.5mm)".

SMI Labelling  Name: <i>Emma Henderson</i> Date: <i>16 May 2014</i> <i>Signature held on master</i>	Regulatory Affairs  Name: <i>A. HOWARD-GRIGG</i> Date: <i>16 May 2014</i> <i>Signature held on master</i>	Marketing  Name: <i>A. HOWARD-GRIGG</i> Date: <i>16 May 2014</i> <i>Signature held on master</i>
Compiled by Dan Johnson 	Design Authority R&D / Operations Technology  Name: <i>A. HOWARD-GRIGG</i> Date: <i>16 May 2014</i> <i>Signature held on master</i>	Date Effective 13 JUN 2014
Compilation Date - 16 MAY 14	Circulation: Czech Republic	Sheet 1 of 19

Folding Diagram

Internal Pages



Back Cover

Front Cover

Internal Flap

Cover Flap Folds inwards


FINISHED FOLDED SIZE A6

 = Stitch locations.
Internal pages stitched twice to cover.

Print Details: PMS Pantone Process Black U.

66

English

en

Portex® Blue Line Ultra®
Tracheostomy Inner Cannula,
Single Patient Use

INSTRUCTIONS FOR USE

These instructions are intended for use with the following Portex® products:
Catalogue No: 100/850/060-100 Replacement smooth inner cannula x 2
Internal REF: 100/850/060CZ-100CZ
Catalogue No: 100/851/060-100 Replacement fenestrated smooth inner cannula x 2
Internal REF: 100/851/060-100CZ

These instructions contain important information for safe use of the product. Read the entire contents of these Instructions For Use, including Warnings and Cautions, before using this product. Failure to properly follow warnings, cautions and instructions could result in death or serious injury to the patient and/or clinician.

NOTE: DISTRIBUTE THIS INSERT SHEET TO ALL PRODUCT LOCATIONS.

Please note that the code number of the Portex® Single Patient Use Inner Cannula indicates the size of Portex® Tracheostomy tube with which it should be used e.g.

Inner Cannula Code No.	For use with Tracheostomy Tube	Inner Cannula I.D.	Jackson Size (approx.)	Tracheostomy Tube Length 'C'
100/850/060 100/851/060	6.0 mm	5.0 mm	5	64.5 mm
100/850/070 100/851/070	7.0 mm	6.0 mm	–	70.0 mm
100/850/075 100/851/075	7.5 mm	6.5 mm	7	73.0 mm
100/850/080 100/851/080	8.0 mm	7.0 mm	8	75.5 mm
100/850/085 100/851/085	8.5 mm	7.5 mm	–	78.0 mm
100/850/090 100/851/090	9.0 mm	8.0 mm	–	81.0 mm
100/850/100 100/851/100	10.0 mm	8.5 mm	10	87.5 mm

NOTE:
1 Use only Blue Line Ultra® Inner Cannula with Blue Line Ultra® Tracheostomy Tubes.
2 The use of an Inner Cannula reduces the airway diameter by 1.0mm (excluding size 10.0mm tracheostomy tube which reduces the airway diameter by 1.5mm).

Description

A range of sterile single patient use replacement inner cannula intended to assist in the provision of a patent airway in the tracheostomised patient, available in both fenestrated (coloured red) and unfenestrated (clear) versions for use with the following Portex® Tracheostomy tubes.

100/800/060-100	100/811/060-100	100/831/070-090	100/860/060-100
100/802/060-100	100/812/060-100	100/832/070-090	100/870/060-100
100/810/060-100	100/813/060-100	100/833/070-090	

Instructions for use

1 Following insertion of the tracheostomy tube and removal of the obturator, insert the inner cannula fully until a click is felt indicating correct positioning.

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fr

Warnings

- Do not lubricate the inner cannula as this may occlude the inner cannula and prevent the inner cannula being retained within the tracheostomy tube.
- Contact with electrosurgery electrodes or laser surgery beams must be avoided because PVC will produce toxic fumes in air, or ignite in an enriched oxygen environment (e.g. anaesthesia).
- Never use a Portex® Fenestrated Tracheostomy tube for manual or mechanical ventilation unless a Portex® non fenestrated inner cannula is in place to prevent gas leakage through the fenestrations in the tracheostomy tube.
- If a Portex® Fenestrated Tracheostomy tube is used with a Decannulation cap, either a fenestrated inner cannula must be used or if a non fenestrated inner cannula is present, it should be removed, and in both cases the tracheostomy tube cuff must be deflated. Failure to do so may lead to partial or complete airway obstruction.
- Caution must be exercised when using a Decannulation Cap in conjunction with a cuffed fenestrated tracheostomy tube. The user must ensure that the cuff is completely deflated and that a fenestrated inner cannula is used, or that the inner cannula has been removed. Failure to do so may result in total or partial obstruction of the airway. Observe the patient for signs of respiratory distress and immediately remove the Decannulation Cap if required.
- Do not lubricate the inner cannula as this may occlude the inner cannula and prevent the inner cannula being retained within the tracheostomy tube.
- If it is not possible to remove the inner cannula from the tracheostomy tube, do not try to forcibly remove. Both the inner cannula and tracheostomy tube should be removed together and replaced with a new tracheostomy tube and inner cannula.
- Ensure that the inner cannula is both the appropriate diameter and length for the tracheostomy tube in use. This can be checked by comparing size markings on the tracheostomy tube flange and patient labels with those on the inner cannula unit pack label. Use this inner cannula only with tracheostomy tubes from the Portex® Blue Line Ultra® range.
- Use of an inner cannula of the wrong diameter may cause difficulties in introduction or unnecessarily restrict gas flow. An inner cannula which is too long may protrude excessively from the outer tube, leading to tracheal damage or occlusion. Use of an inner cannula which is too short may lead to build-up of secretions which could cause infection or blockage.
- The use of an Inner Cannula reduces the airway diameter by 1.0mm (excluding size 10.0mm tracheostomy tube which reduces the airway diameter by 1.5mm).
- Do not fit an inner cannula when using a tracheostomy tube in conjunction with a right-angled connector that incorporates a co-axial fresh gas tube which protrudes into the connector lumen (e.g. Norman elbow). The protruding fresh gas tube could occlude the inner cannula, restricting expiratory flow of gas, resulting in barotrauma/hypoxia.
- During and after attachment of the breathing system to the tracheostomy tube connector, avoid application of excessive rotational or linear forces on the tube to prevent accidental disconnection or occlusion.

Portex, Blue Line Ultra, Thermavent, Portex design mark and Smiths Medical design mark are trademarks of Smiths Medical. The symbol ® indicates the trademark is registered in the U.S. Patent and Trademark office and certain other countries. All other names and marks mentioned are the trade names, trademarks or service marks of their respective owners.

© 2014 Smiths Medical. All rights reserved.

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Portex, Blue Line Ultra, Thermovent designmerke, Portex designmerke, samt Smiths Medical designmerke er varemerker som tilhører Smiths Medical. Symbolet ® angir at varemerket er registrert ved U.S. Patent and Trademark Office og i visse andre land. Alle andre navn og merker som er nevnt, er varemerker, handelsnavn og servicemerker som tilhører de respektive elere.