**Cutimed® Siltec® Sorbact® B /  
Cutimed® Siltec® Sorbact® B Sacrum**

Read the instructions for use before starting treatment with Cutimed® Siltec® Sorbact® B / Cutimed® Siltec® Sorbact® B Sacrum.

**Device Description**

Cutimed® Siltec® Sorbact® B / Cutimed® Siltec® Sorbact® B Sacrum is a hydrophobic microbe binding wound dressing, based on Sorbact® Technology. It consists of a Sorbact® DACCTM-coated, wound contact layer combined with an absorbent polyurethane foam that contains super-absorbent stripes and a border with silicone adhesive, which minimises trauma during dressing changes. Cutimed® Siltec® Sorbact® B / Cutimed® Siltec® Sorbact® B Sacrum manages wound infections. The outer film is water-repellent, yet permeable to oxygen and vapour. The foam dressing and the super-absorbent stripes are designed to absorb and to lock away excess wound fluid inside the dressing, promoting a moist wound environment and minimizing the risk of maceration. Cutimed® Siltec® Sorbact® B / Cutimed® Siltec® Sorbact® B Sacrum is highly conformable and provides protection of the wound site. Cutimed® Siltec® Sorbact® B / Cutimed® Siltec® Sorbact® B Sacrum leaves no residues in the wound. The dressing reliably retains wound exudate, even under compression.

**Mode of Action**

Sorbact® can bind hydrophobic microbes based on *in vitro* testing of the following species: *S. aureus* including MRSA, *S. gordonii*, *P. aeruginosa*, and *E. coli*. When the dressing is removed, microbes adhered to the dressing will be removed.

**Intended Use**

Cutimed® Siltec® Sorbact® B / Cutimed® Siltec® Sorbact® B Sacrum is indicated for use in the management of low to high exuding partial to full thickness wounds (including clean, colonized, contaminated, or infected wounds): chronic wounds (venous and arterial ulcers, diabetic ulcers and pressure ulcers), postoperative dehisced wounds and traumatic wounds.

**Contraindications**

There are no known contraindications to the use of Cutimed® Siltec® Sorbact® B / Cutimed® Siltec® Sorbact® B Sacrum. Warnings and precautions should be noted.

**Warnings and Precautions**

Do not re-use. Cutimed® Siltec® Sorbact® B / Cutimed® Siltec® Sorbact® B Sacrum is for single patient and single use only. Re-use may lead to cross contamination and increased risk of infection. Do not use if the peel pouch is already opened or damaged, since the sterility cannot then be guaranteed. Do not use on patients with a known hypersensitivity to the dressing (dialkylcarbamoyl chloride [DACCTM], cellulose acetate, polyurethane, silicone and sodium polyacrylate). Do not use in combination with fatty products, such as ointments, creams and solutions, as they may decrease the binding of hydrophobic microbes. Do not use with oxidizing solutions such as hypochlorite or hydrogen peroxide. Do not re-sterilize. If the treated condition deteriorates, fails to improve or a side effect is observed, consult an appropriate clinician.

**Instructions for Use**


1. Prepare the wound and surrounding skin according to local clinical practice. Ensure that the surrounding skin is clean and dry.
  2. Select an appropriate dressing size and shape that covers the wound and is suitable for the anatomical location. The wound pad should overlap the wound margins. (A) Do not cut the dressing.
  3. Remove the dressing from the peel pouch using an aseptic technique.
  4. Remove the protective film (B) from the dressing and apply the dressing to the wound. Ensure that the green Sorbact® wound contact layer comes into direct contact with the complete wound surface to allow hydrophobic microbes to bind to the dressing. (C)
  5. Avoid stretching the dressing.
  6. The dressing change frequency depends on exudate levels and the overall condition of the wound and surrounding skin. Should the clinical condition allow, the dressing can be left in place for up to 7 days.
  7. Gently remove Cutimed® Siltec® Sorbact® B / Cutimed® Siltec® Sorbact® B Sacrum from the wound and discard.
- Cutimed® Siltec® Sorbact® B / Cutimed® Siltec® Sorbact® B Sacrum can be used in conjunction with compression therapy.
- Cutimed® Siltec® Sorbact® B / Cutimed® Siltec® Sorbact® B Sacrum is MRI safe.
- Prior to radiation therapy, remove Cutimed® Siltec® Sorbact® B / Cutimed® Siltec® Sorbact® B Sacrum if placed in the radiation field. A new dressing can be applied following treatment.

**Storage and Disposal**

Cutimed® Siltec® Sorbact® B / Cutimed® Siltec® Sorbact® B Sacrum shall be stored dry and kept away from sunlight. Disposal should be made according to local environmental procedures.

**Caution**

Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

**Rx only**

|                                |                                                      |                                   |                         |                                                                       |
|--------------------------------|------------------------------------------------------|-----------------------------------|-------------------------|-----------------------------------------------------------------------|
| <b>REF</b><br>Catalogue number | <b>LOT</b><br>Batch code                             | <br>Use-by date                   | <br>Date of manufacture | <b>UDI</b><br>Unique device identifier                                |
| <b>MD</b><br>Medical device    | <br>Caution                                          | <br>Consult instructions for use  | <br>Do not re-use       | <br>Keep away from sunlight                                           |
| <br>Keep dry                   | <b>STERILE EO</b><br>Sterilized using ethylene oxide | <br>Single sterile barrier system | <br>Do not re-sterilize | <br>Do not use if package is damaged and consult instructions for use |

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Cutimed® is a registered trademark of BSN medical GmbH.
Siltec® is a registered trademark of BSN medical GmbH.
Sorbact® is a registered trademark of ABIGO Medical AB.
Cutimed® es una marca comercial registrada de BSN medical GmbH.
Siltec® es una marca comercial registrada de BSN medical GmbH.
Sorbact® es una marca comercial registrada de ABIGO Medical AB.