[illegible]

Para información sobre ALEVEEN LIFE
para más informaciones sobre el
producto ir al código QR por maggot
NLIFE / Scan on QR-code voor meer
informatie over ALEVEEN LIFE /
Informazioni su ALEVEEN LIFE
Information on ALEVEEN LIFE /
ALEVEEN LIFE / Oskovinarit QR koodi
Täza dodatne informácie o výrobku
QR-koodi / Запознате твој код QR
IF / Z ALIVEEEN LIFE je varnostno
prijemljiva za ALEVEEN LIFE, naskalnikite
na ALEVEEN LIFE s nalepkami razlozavate ved QR
aleveen life / Naskenken QR-
code informatie o ALEVEEN LIFE oddajte
lokator.

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Instructions for Use: Select the appropriate dressing size, ensure
the entire wound or area to be protected. The scale on the skin
only! Cleanse the application site. Ensure surrounding skin is clean
and free of excess hair.

1.Remove first liner and
anchor the adhesive side of
the dressing to the skin.

2.Place the narrow end at
a minimum of 2cm (3/4")
above the anal sphincter,
then smooth the dressing
over the sacrum.

3.Cleanse the skin around the
dressing edge.
4.Apply additional dressing
if needed.
5.Repeat steps 1-4 as
needed.

ALEVEEN LIFE Dressings can be used in conjunction with INTRA
neutrotic or sloughy wounds.

Observation: When using this pressure injury prevention, frequent
pressure damage by lifting back one corner, inspecting and then
Removal: Lift the top edge and slowly peel downwards towards the
bottom. Do not pull straight down. This will minimise the chance of trauma
to the skin. Improper disposal of dressing and packaging can cause slips, cuts
or infection; dispose of them in accordance with local clinical protocols or hygiene
guidelines.

Frequency of change: During the early stages of treatment, inspect
frequently. Sacral dressings can be left in place for up to 7 days,
changed depending on the condition. Change if there are signs of
change: 75% of the change indicator. Consider changing if the ex-
cess of 50% of the change indicator. The diagram below is for guidance
increased potential for contamination and infection in the presence of
dressing adherence may be required as per local clinical protocol
healthcare professional.

1. No need to change 2. No need to change 3. Consider change

© Smith+Nephew Ltd. All rights reserved. 2025

Trade Marks of Smith+Nephew, Date of issue 11/2025, Revision 01
Covered by one or more patents (www.smith-nephew.com)

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