

PRODUCT DATASHEET

Disposable Slipper (Fracture) Pan Liner


1.7L, Medium (Midi), w/Lid, WardPak™

Brand	: CURAS®
Category	: Nursing
Range	: Hygiene Pulp™
Product Group	: Human Waste Collection
REF	: 08001



PRODUCT DESCRIPTION

Disposable slipper (fracture)pan liner made of pulp material with water resistance properties. 1.7L capacity. Suitable for disposal in a macerator. Lids are included for emergency use. Packed in WardPak™ for easy dispensing. Non-sterile.



INTENDED USE

Disposable slipper (fracture)pan liner for holding human waste and other body fluids. Also used by immobile female patients for bed-side collection of urine. Can be used in conjunction with a standard plastic support for improved stability. Single use only.

FEATURES& BENEFITS

- An environmental friendly product that is made of 100% recycled newsprint paper.
- Disposable for reduced risk of cross-contamination and Hospital-Acquired Infections.
- Water resistant, for a minimum period of 4 hours (applicable to slipper pan only).
- Integrated graduation marking for indication of approximate volume.
- Stack 2 pieces or more for extra stability and improved hygiene, if required.
- Fits perfectly into standard plastic support (REF 08111) for extra stability.
- Disposable lid helps to contain odours and prevent accidental spills.
- Stackable for reduced storage space.
- Steady with flat bottom.
- Inner bag, WardPak™ protects product from dust and contamination.

PRODUCT SPECIFICATIONS

REF	Capacity	Dimensions (L x W x H)	Average Weight/Pc	Packing
08001	1.7L	34.3 x 29.6x8.4 cm	56g ± 4 g	WardPak™

Composition

- Slipper (fracture) pan liner is made of 100% recycled newsprint paper with water-based fluid resistant compound.
- Lid is made of cardboard.

Country of Origin

- China

Sterilization

- Non-sterile.
- Can be sterilized by EtO and Gamma radiation.

Shelf Life

- Slipper (Fracture) Pan Liner & Lid: 8 years

Selected Standards & Guidelines

- MDD 93/42/EEC (CE Marking, Class I), as amended by Directive 2007/47/EC.
- EN ISO 9001:2008
- PAS29:1999 Fluid resistance test for a minimum period of 4 hours (water at 35°C±3°C).
- PAS29:1999 Maceration test with a hospital sluice room type macerator, normal cycle.

Recommended Storage Condition

- Store at room temperature and avoid exposure to extreme temperatures.
- Keep at cool dry place, away from moisture, heat and sources of ignition.

Proposed Disposal Method

- If regular waste, dispose of used product into a macerator, or in accordance with the local waste protocol. If a macerator is not available, incinerate or land-fill.
- If content is contaminated or contagious, treat as clinical waste and dispose of in accordance with the clinical waste protocol.

Product Symbols



PACKAGING INFO

REF	EAN-13 Barcode	WardPak™	Carton	Gross Weight
08001	5051400080012	50 pcs / bag	Slipper (Fracture) Pan Liner: 150 pcs / carton Lid : 5 pcs / carton	~9 kg

Components / Material Properties			Recommendations on Waste Handling				
Packaging	Material	Contains PVC	Incinerate	Landfill	Macerate	Recycle	Reuse
WardPak™	PE	--	√	(√)		√	
Outer Carton	Corrugated Carton	--	(√)	(√)		√	
Transport Pallet	Wood	--				(√)	√

√ = Highly Recommended (√) = Recommended

ADDITIONAL INFO

- Curas® Hygiene Pulp™ products have been tested in several macerator models and are compliant under recommended operating procedures.
- HS Code: 4823709000
- Documentation available upon request. Kindly contact: info@curas.com