

PRODUCT DATASHEET

Disposable Urine Bottle, Round 1L, w/SAP by Curas™ & Lid, RoomPak™

Brand	: CURAS®
Category	: Nursing
Range	: Hygiene Pulp™
Product Group	: Disposable Urine Bottles
REF	: 08058



PRODUCT DESCRIPTION

Disposable round urine bottle with lid made of pulp material with water resistance properties. 1L capacity. Suitable for disposal in a macerator. Non-sterile. Packed in RoomPak™, a type of dispensing bag for easy dispensing and enhance hygienic performance. Product comes with Curas™ SAP (Super Absorbent Polymer) powder to solidify urine into gel to prevent spilling and reduce odour as well as improve waste handling.

INTENDED USE

Disposable round urine bottle for use by immobile male patients for collection of urine. The bottle is also ideal for collecting urine drained from urine bag. Suitable for single use.

FEATURES & BENEFITS

- Disposable and single use can reduce risk of cross-contamination and Hospital-Acquired Infections (HAI).
- Water resistant, for a minimum period of 4 hours (applicable to urine bottle and lid only).
- Ergonomic jug design for convenient use and holding optimal volume.
- Reinforced flat base provides extra strength and allows stable placement.
- Designed to easily fit into bedside urine bottle holder.
- Disposable lid helps to control odours and accidental spills.
- All in one package for convenient one-stop access to the necessary related products.
- RoomPak™ dispensing bag protects urine bottle from dust and contamination.
- Perforated opening on dispensing bag offers an easy dispensing option.
- Comes with super absorbent polymer (SAP) by Curas™ to prevent accidental spills, control odour and improve waste handling.

PRODUCT SPECIFICATIONS

REF	Capacity	Dimensions (L x W x H)	Average Weight/Pc	Packing
08058	1L	24.5 x 11 x 12 ± cm	Urine bottle = 50 g ± 5 g Lid = 3-4 g SAP = ~9 g	RoomPak™

Composition

- Urine bottle and lid made of 100% recycled newsprint paper with water-based fluid resistant compound.
- Superabsorbent copolymer (sodium polyacrylate) powder.
- Physical state: White granular solid that is odourless.
- Water-soluble sachet (pouch).

Country of Origin

- China

Sterilization

- Non-Sterile
- Can be sterilized by EtO and Gamma radiation.

Shelf Life

- Urine bottle = 8 years
- SAP = 5 years

Selected Standards & Guidelines

- MDD 93/42/EEC (CE Marking, Class I), as amended by Directive 2007/47/EC.
- ISO 9001:2008
- ISO 2859-1 Sampling Plan
- PAS29:1999 Fluid resistance test for a minimum period of 4 hours (water at $35^{\circ}\text{C} \pm 3^{\circ}\text{C}$).
- PAS29:1999 Maceration test with a hospital sluice room type macerator, normal cycle.

Recommended Storage Condition

- Store at room temperature and avoid exposure to extreme temperatures.
- Keep at cool dry place, away from moisture, heat and sources of ignition.

Proposed Disposal Method

- If regular waste, dispose of used product together with content into a pulp macerator, or in accordance with the local waste protocol. If a macerator is not available, incinerate or land-fill.
- If content is contaminated or contagious, treat as clinical waste and dispose of in accordance with the clinical waste protocol

Product Symbols

curas.com/Infection-Prevention



PACKAGING INFO

REF	EAN-13 Barcode	RoomPak™	Carton	Gross Weight
08058	5051400080586	Urine Bottle : 5 pcs / bag	Urine Bottle : 40 pcs / carton	~3.6 kg
		Lids: 5 pcs / bag	Lids : 40 pcs / carton	
		SAP: 5 pcs / bag	SAP : 40 pcs / carton	

Components / Material Properties			Recommendations on Waste Handling				
Packaging	Material	Contains PVC	Incinerate	Landfill	Macerate	Recycle	Reuse
RoomPak™	PE bag	--	√	(√)		√	
Carton	Brown Corrugated Carton	--	(√)	(√)		√	
Transport Pallet	Wood	--				(√)	√

√ = Highly recommended

(√) = Recommended

ADDITIONAL INFO

- Curas® Hygiene Pulp™ products have been tested in several macerator models and are compliant under recommended operating procedures.
- HS Code: 4823709000
- Documentation available upon request. Kindly contact: info@curas.com