

PRODUCT DATASHEET

Disposable Urine Bottle, Round 1L, w/Lid & Measuring Stick

Brand	: CURAS®
Category	: Nursing
Range	: Hygiene Pulp™
Product Group	: Disposable Urine Bottles
REF	: 08050



PRODUCT DESCRIPTION

Disposable round urine bottle made of hygienic pulp material with water resistance properties. 1L capacity. Comes with complimentary lids and measuring sticks. Suitable for disposal in a macerator.



INTENDED USE

Disposable round urine bottle for use by immobile male patients for collection of urine and draining of urine bags. Single use.

FEATURES & BENEFITS

- Disposable for reduced risk of cross-contamination and Hospital-Acquired Infections.
- Disposable lid helps to contain odours and prevent accidental spills.
- Disposable measuring stick helps to measure volume of urine.
- Product is steady with flat bottom.
- Fits into bedside urine bottle holder.
- Water resistant, for a minimum period of 4 hours (applicable to urine bottle and lid only).
- Fluid Absorber, SAP (REF 08100) can be inserted prior to use for improved waste handling.
- Dispenser-style carton held off the floor and away from contamination and infections.

PRODUCT SPECIFICATIONS

REF	Capacity	Dimensions (L x W x H)	Average Weight/Pc	Packing
08050	1L	24.5 x 11 x 12 cm	Urine bottle: 50 ± 5 g Lid: 3 – 4 g Stick: ~2.5 g	Dispenser Carton

Composition

- Urine bottle and lid made of 100% recycled newsprint paper and water-based fluid resistant compound.
- Measuring stick made of cardboard.

Country of Origin

- China

Sterilization

- Non-Sterile
- Can be sterilized by EtO and Gamma radiation.

Shelf Life

- 8 years

Selected Standards & Guidelines

- MDD 93/42/EEC (CE Marking, Class I), as amended by Directive 2007/47/EC.
- ISO 9001:2008
- PAS29:1999 Fluid resistance test for a minimum period of 4 hours (water at 35°C±3°C).
- PAS29:1999 Maceration test with a hospital sluice room type macerator, normal cycle.

Recommended Storage Condition

- Store at room temperature and avoid exposure to extreme temperatures.
- Keep at cool dry place, away from moisture, heat and sources of ignition.

Proposed Disposal Method

- If regular waste, dispose of used product together with lid and measuring stick into a macerator, or in accordance with the local waste protocol. If a macerator is not available, incinerate or land-fill.
- If content is contaminated or contagious, treat as clinical waste and dispose of in accordance with the clinical waste protocol.

Product Symbols



curas.com/Symbols

PACKAGING INFO

REF	EAN-13 Barcode	Dispenser Carton	Gross Weight
08050	5051400080500	Urine Bottles: 120 pcs / carton Lids: 10 pcs / carton Measuring Sticks: 5 pcs / carton	~7.9 kg

Components / Material Properties			Recommendations on Waste Handling				
Packaging	Material	Contains PVC	Incinerate	Landfill	Macerate	Recycle	Reuse
Dispenser Carton	Corrugated Carton	--	(√)	(√)		√	
Transport Pallet	Wood	--				(√)	√

√ = Highly recommended

(√) = Recommended

ADDITIONAL INFO

- Curas® Hygiene Pulp™ products have been tested in several macerator models and are compliant under recommended operating procedures.
- HS Code: 4823709000
- Documentation available upon request. Kindly contact: info@curas.com