

PRODUCT DATASHEET

Disposal Deep Trays

Brand	: CURAS [®]
Category	: Nursing
Range	: Hygiene Pulp [™]
Product REF	: 08202, 08203, 08204
Product Description	: Disposable Deep Tray, 100% Pure Natural Fibre, WardPak [™]



PRODUCT DESCRIPTION

Disposable deep trays made of pure natural fibre with water, heat, steam and oil resistant properties. Biodegradable, compostable and suitable for disposal in a macerator*. Packed in WardPak[™] (inner bag) for easy dispensing at ward level. Non-sterile.

INTENDED USE

Fibre trays are used for various applications, specifically at the ward level, isolation ward, ICU and in the operating room. They are used for holding and transporting drugs and devices such as swabs, gauze, syringes, dressing, syringes, needles, instruments, and as trays in procedure kits. They can also be used to hold wipes and other products for heating in microwave or oven. These trays are ideal for sterile procedure kit packs, since they can be autoclaved.

The deeper trays (size M and above) can be used to hold waste arising from medical, nursing or dental procedures, e.g., human tissue, blood or other body fluids. Single use only.

FEATURES & BENEFITS

- An environmental friendly product that is biodegradable, and compostable.
- Clean, hygienic white trays made of 100% natural fibre.
- Disposable and single use can reduce risk of cross-contamination and Hospital-Acquired Infections.
- ECF (elemental chlorine free) bleaching method for a safe colour free product.
- Dust-free and non-toxic; safe for medical use and food contact.
- UV-treated for improved surface hygiene and cleanliness.
- Safe for use in autoclave machine, microwave, oven and refrigerator.
- Autoclavable for procedure kit preparation (120°C to 136°C)
- Resistant to hot water (up to 100°C), hot oil (up to 120°C), and steam (up to 136°C).
- No toxic substance or odour is released in high-temperature or acidic/alkaline conditions.
- Smooth, stackable and easy to separate.
- Packed in WardPak[™] (inner bag) to protect products from dust and contamination.

PRODUCT SPECIFICATIONS

REF	Size	Dimensions (L x W x H)	Average Weight/Pc	WardPak [™]
08202	S	18x12.5x3.5cm(7.1x4.9x1.4in)	10 g	Inner Bag
08203	M	20.5 x12.5 x 4.5cm(8.1x4.9x1.8in)	16 g	Inner Bag
08204	L	27x19x4.5cm (10.6x7.5x1.8in)	30 g	Inner Bag

Composition

- Made of 100% natural fibre material (sugarcane bagasse).
- ECF (elemental chlorine free) bleaching method.
- Food-grade fluid and oil-resistant compound.

Country of Origin

- China

Sterilization

- Non-sterile.
- Can be sterilized by EtO, Gamma radiation and autoclave.

Shelf Life

- 8 years

Selected Standards & Guidelines

- MDD93/42/EEC (CE Marking, Class I), as amended by Directive 2007/47/EC
- ISO 9001:2008
- ISO 14001
- EN 13432:2000/12
- ASTM D-6868
- US FDA 21 CFR 176.170
- BRC Global Standards
- HACCP
- PAS29:1999 Fluid Resistance test for a minimum of 4 hours (water at 35°C±3°C)
- PAS29:1999 Maceration test*

Recommended Storage Condition

- Store at room temperature and avoid exposure to extreme temperatures.
- Keep at cool dry place, away from heat and moisture.

Proposed Disposal Method

- If regular waste, dispose of used product into a macerator, or in accordance with the local waste protocol. If a macerator is not available, incinerate or land-fill.
- If content is contaminated or contagious, treat as clinical waste and dispose of in accordance with the clinical waste protocol.

Product Symbols



PACKAGING INFO

REF	EAN-13 Barcode	WardPak™	Carton	Gross Weight
08202	5051400082023	100 pcs / bag	1,200 pcs / carton	~13 kg
08203	5051400082030	100 pcs / bag	800 pcs / carton	~14 kg
08204	5051400082047	50 pcs / bag	500 pcs / carton	~16 kg

Components / Material Properties			Recommendations on Waste Handling				
Packaging	Material	Contains PVC	Incinerate	Landfill	Macerate	Recycle	Reuse
WardPak™	PE	--	√	(√)		√	
Outer Carton	Corrugated Carton	--	(√)	(√)		√	
Transport Pallet	Wood	--				(√)	√

√ = Highly recommended

(√) = Recommended

ADDITIONAL INFO

- Curas® Hygiene Pulp™ products have been tested in several macerator*models and are compliant under recommended operating procedures.
- HS Code: 4823709000
- Documentationavailable upon request. Kindly contact: info@curas.com

** Some macerators will not accept these Trays. Others may require more than one cycle. Consult your local technician.*