

PRODUCT DATASHEET

Disposable Urine Bottle, Round 1L, w/SAP by Curas™, Lid & Stick, WardPak™

Brand	: CURAS®
Category	: Nursing
Range	: Hygiene Pulp™
Product Group	: Disposable Urine Bottles
REF	: 08056



PRODUCT DESCRIPTION

Disposable round urine bottle with SAP by Curas™ and Lid. Comes with complimentary graduated sticks for approximate measurement. Suitable for disposal in a macerator. 1L capacity.

INTENDED USE

Disposable round urine bottle for use by immobile male patients for collection of urine and draining of urine bags. Packed with SAP by Curas™ and lid to prevent spills, improve waste handling, and control odour. Single use.



FEATURES & BENEFITS

- Disposable and single use can reduce risk of cross-contamination and Hospital-Acquired Infections.
- Urine bottle is designed with a flat bottom for stability and easy placement.
- Ergonomic jug design for convenient use and holding optimal volume.
- Able to fit into common bedside holders.
- Water resistant, for a minimum period of 4 hours (applicable to urine bottle and lid only).
- Preloaded with SAP by Curas™ to crystalize liquid on contact for convenience and safe handling.
- SAP by Curas™ and disposable lid helps to contain odours and prevent accidental spills.
- Disposable graduated stick provides approximate measure of urine collected.
- All in one package for convenient one-stop access to the necessary related products.
- WardPak™, dispensing carton protects product from dust and contamination.
- Perforated opening on dispensing carton offers an easy dispensing option.

PRODUCT SPECIFICATIONS

REF	Capacity	Dimensions (L x W x H)	Average Weight/Pc	Packing
08056	1L	24.5 x 11 x 12 cm	Urine bottle: 50 ± 5g Lid: 3 – 4 g Stick: ~2.5 g SAP: ~9 g Bottle + Lid + Stick + SAP: 64.5 g	WardPak™

Composition

- Urine bottle and lid made of 100% recycled newsprint paper with water-based fluid resistant compound.
- SAP by Curas™ made of superabsorbent copolymer (sodium polyacrylate) powder.
- Graduated stick made of cardboard.

Country of Origin

- China

Sterilization

- Non-Sterile
- Can be sterilized by EtO and Gamma radiation.

Shelf Life

- 5 years for SAP.
- 8 years for urine bottle, lid and graduated stick.

Selected Standards & Guidelines

- MDD 93/42/EEC (CE Marking, Class I), as amended by Directive 2007/47/EC.
- ISO 9001:2008
- PAS29:1999 Fluid resistance test for a minimum period of 4 hours (water at 35°C±3°C).
- PAS29:1999 Maceration test with a hospital sluice room type macerator, normal cycle.

Recommended Storage Condition

- Store at room temperature and avoid exposure to extreme temperatures.
- Keep at cool dry place, away from moisture, heat and sources of ignition.

Proposed Disposal Method

- If regular waste, dispose of used product together with content and SAP (solidified liquid) into a macerator, or in accordance with the local waste protocol. If a macerator is not available, incinerate or land-fill.
- If content is contaminated or contagious, treat as clinical waste and dispose of in accordance with the clinical waste protocol.

Product Symbols

Urine Bottle Symbols



curas.com/Symbols

SAP by Curas™

Symbols



PACKAGING INFO

REF	EAN-13 Barcode	WardPak™ (pcs / pack)	Gross Weight
08056	5051400080562	Urine Bottles : 28pcs / carton SAP : 28pcs / carton Lids : 28 pcs / carton Graduated Sticks : 10 pcs / carton	~2.3 kg

Components / Material Properties			Recommendations on Waste Handling				
Packaging	Material	Contains PVC	Incinerate	Landfill	Macerate	Recycle	Reuse
WardPak™	White Corrugated Carton	--	(√)	(√)		√	
Transport Pallet	Wood	--				(√)	√

√ = Highly recommended

(√) = Recommended

ADDITIONAL INFO

- Curas® Hygiene Pulp™ products have been tested in several macerator models and are compliant under recommended operating procedures.
- HS Code: 4823709000
- Documentation available upon request. Kindly contact: info@curas.com