

carefix[®]

Creating functional comfort



FINGER & TOE fixation

Precut elastic dressing fixation



TYTEX[®]

The CareFix finger/toe fixation is designed in soft and elastic material, that makes it easy to use and comfortable to the patient. It can be applied to the finger as well as the toe.

Features:

- One product for both finger and toe fixation
- Ready to use. No tools of any kind needed
- Simple and quick application
- High sustained elasticity
- Fixes dressings and bandages
- Comfortable for the patient
- Practical colour-coding indicates sizes
- Non-irritant to the skin
- Washable
- Latex-free material for the sake of the patient and caregiver
- Tested for harmful substances according to Oeko-Tex standard 100

Usage after treatment of injury on finger or toe

The fixation can be worn alone or can be used to hold surgical dressings in place without tape.

The CareFix finger/toe provides extra security to the wound. It reduces the risk of dirt in the wound, which again reduces the risk of infection, while reminding the patient of the injured finger/toe.

Composition:

79% Polyamide, 21% Elastane.

Washing instructions:



60°C/140°F


Recommended sizes:

S = Child
M = Adult
L = Adult

Packaging:



Dispenser boxes w/20 pcs.



Application

For further information please contact:

Dealer:

DS-Finger&Toe-01-GB-0515